# Universal open platform and reference specification for Ambient Assisted Living



# Personal Health Devices (x073) and universAAL middleware integration

<< How to for uAAL developers >>

Local device discovery and integration group (LDDI)

## **Release history**

Version	Date	Editors Reviewers		Comments
0	July 31, 2012	Ángel Martínez (TSB) and Luis Gigante (TSB)	Juan-Pablo Lázaro (TSB)	1st version
1	September 17, 2012	Ángel Martínez (TSB)	Thomas Fuxreiter (AIT) and Patrick Stern (AIT)	2nd version
2	December 10, 2012	Ángel Martínez (TSB)	Thomas Fuxreiter (AIT) and Alejandro Medrano (UPM)	3rd version

### **Table of contents**

3
6
6
8
8
8
10
16
17
21
29
31
31
39

#### **Overview**

The aim of this document is to make available a technical guide to enable universAAL (uAAL) developers community in general and LDDI experts group in particular to start developing Ambient Assisted Living (AAL) services with Continua Health Alliance devices (also called agents or sources).

This technological solution involves two pieces of code:

- libHDPnative.so is a dynamic library fully developed in C code that enables communication between BlueZ Bluetooth stack (C) and uAAL services (Java) through Java Native Interface (JNI) technology and interprocess communication (IPC) DBUS system. Is an isolated file with ".so" extension (similar to DLL files in Microsoft Windows environments).
- x73implementation.jar is the Continua Health Alliance standard (ISO/IEEE 11073) implementation developed in pure Java and without any external dependencies. It is composed by several packages grouped by performance: coder, decoder, DIM, state machine, etc.

Additionally other files or packages will be distributed:

- javadoc.rar file with a short description of methods, attributes, classes and relationships between them. Quite simply for Java developers.
- A sample project (BloodPressureMonitor\_x073\_Sample) showing how a
  Continua manager can receive valid and real time measurements
  through a HDP wireless channel. This example works outside OSGi
  platforms (no uAAL middleware is needed) so it should just be consider
  as a proof of concept.
- Some Maven projects showing how a Continua manager can receive weighting scales as well as blood pressure monitor events and publish it inside uAAL context bus is provided. Both sample applications can be followed as a template in order to develop additional managers with more advance features.

Just to put in context the future reader, this guide (text and images) have been elaborated with the following test environment (not relevant at this point):

- Sony Vaio laptop with Ubuntu 11.10 OS (kernel version 3.0.0-19)¹.
- GNU/Linux bluetooth stack blueZ<sup>2</sup> 4.99 version.
- Anycom 2.0 USB-200 bluetooth dongle.
- Real Continua health devices (just for testing): UC-321PBT-C A&D weigh scale<sup>3</sup> and UA-767PBT-C A &D upper arm blood pressure monitor<sup>4</sup>.

<sup>1</sup> OS had successfully been updated to ubuntu 12.04 LTS in August 2012.

<sup>2 &</sup>lt;a href="http://www.bluez.org">http://www.bluez.org</a>

<sup>3 &</sup>lt;a href="http://www.aandd.jp/products/medical/bluetooth/uc\_321pbt\_c.html">http://www.aandd.jp/products/medical/bluetooth/uc\_321pbt\_c.html</a>

<sup>4 &</sup>lt;a href="http://www.aandd.jp/products/medical/bluetooth/ua\_767pbt\_c.html">http://www.aandd.jp/products/medical/bluetooth/ua\_767pbt\_c.html</a>

#### **Execution environment**

As has previously been mentioned, this technological solution relies over the official bluetooth stack (BlueZ) available in GNU/Linux operating systems, so from this point you should continue working with any free and open source OS distribution. Due to its high population and success, we decide work with the **version 11.10 of Ubuntu Desktop OS (32 bits)** directly downloaded from Ubuntu web page<sup>5</sup> (Ubuntu 12.04 LTS is out now).



First step is ensure that your computer has the minimum technological requirements necessary to run and start developing new AAL services with the aforementioned libraries:

#### **Hardware requirements**

As Continua health devices employs bluetooth as transmission system basis, in order to send and receive x073 data frames between agents and managers it is obvious that your PC should support that wireless technology. **Developed source code is independent of any bluetooth dongle, so it doesn't matter which model or manufacturer you choose to test the libraries**.

In our case, we choose a bluetooth dongle USB-200 model (Bluetooth 2.0 + Enhanced Data Rate (EDR)) manufactured in 2006 by ANYCOM just because it was available for us as free resource in TSB company (please remember: any technical or strategic reason justifies this dongle choice).



<sup>5 &</sup>lt;a href="http://www.ubuntu.com/download/desktop">http://www.ubuntu.com/download/desktop</a>

In order to verify that your bluetooth dongle works fine in GNU/Linux platforms you can use (for instance) sdptool (interface for performing service queries on bluetooth devices around us) command as follows:

```
$ sdptool browse
Inquiring ...
Browsing 4C:ED:DE:0B:0D:98 ...
Service Name: Audio Source
Service Provider: TOSHIBA CORPORATION
Service RecHandle: 0x1008e
Service Class ID List:
 "Audio Source" (0x110a)
Protocol Descriptor List:
  "L2CAP" (0x0100)
   PSM: 25
  "AVDTP" (0x0019)
   uint16: 0x102
Language Base Attr List:
 code ISO639: 0x656e
  encoding: 0x6a
 base offset: 0x100
Profile Descriptor List:
  "Advanced Audio" (0x110d)
    Version: 0x0102 ...
```

Table 1: Ubuntu console (uAAL developer input/system output)

```
angel@angel-laptop:~$ sdptool browse
Inquiring ...
Browsing 4C:ED:DE:0B:0D:98 ...
Service Name: Audio Source
Service Provider: TOSHIBA CORPORATION
Service RecHandle: 0x1008e
Service Class ID List:
  "Audio Source" (0x110a)
Protocol Descriptor List:
  "L2CAP" (0x0100)
    PSM: 25
  "AVDTP" (0x0019)
    uint16: 0x102
Language Base Attr List:
 code_IS0639: 0x656e
 encoding:
              0хба
  base offset: 0x100
Profile Descriptor List:
  "Advanced Audio" (0x110d)
   Version: 0x0102
```

Figure 1: Available bluetooth services in our coverage range (checking bluetooth dongle performance)

Important note: sdptool command is part of BlueZ (bluetooth tools and daemons) package so if that package is not available in our PC, we can't use it. Just ensure that bluez is ready to be launched or install it from any software manager (Synaptic, Ubuntu software manager centre, etc) otherwise.




Figure 2: Ubuntu software manager with BlueZ package installed

#### **Software requirements**

This technical guide assumes that you have some expertise developing uAAL services as well as setting up the required environment and related plug-ins (if not, please carefully read the universAAL Developer Handbook<sup>6</sup> (login required) and come back later), so you will probably have within your reach some generic Java developer tools, namely Maven, SVN client, Pax Runner and Eclipse IDE. Sorry for bothering you again in that case:

• The Java JDK available by default in GNU/Linux systems is a FLOSS (free and open source) alternative named as OpenJDK. This libraries has been developed with one of the latest versions of the official Oracle JDK available for download so it is highly recommended work with the same JDK than us. After completely remove the OpenJDK from our system (if necessary) we can download and install the right version of the Java Development Kit. Next table shows you how install Oracle-jdk7 for 32 bit architecture (Important note: 64 bit Oracle JDK version don't work on 32

bit OS, ensure you are working with the appropriate binaries package).

```
$ sudo add-apt-repository ppa:webupd8team/java
$ sudo apt-get update
$ sudo apt-get install oracle-jdk7-installer
$ java -version
java version "1.7.0_04"
Java(TM) SE Runtime Environment (build 1.7.0_04-b20)
Java HotSpot(TM) Server VM (build 23.0-b21, mixed mode)
```

Table 2: Downloading and installing Oracle JDK version 7 (32 bits)

Eclipse IDE and related plugins suggested to improve and facilitate uAAL development working process. Please, read <u>Appendix A</u> at the end of this document to a better understanding of how to proceed at each case. Summarizing, you will download **Eclipse IDE for GNU/Linux OS version 3.5.2** (as recommended in the Developer Handbook of universAAL) and install some additional software tools.



Figure 3: Eclipse IDE installed

- As has previously been mentioned, BlueZ is the core of this technological solution and contains tools and system daemons for using Bluetooth devices. BlueZ is the official Linux Bluetooth protocol stack, which has full support of Health Device Profile (HDP) required to exchange data between Continua agents and managers (AAL software services running inside OSGi framework, for instance). Despite of this code is distributed by default in GNU/Linux OS, as developers we need to install next package in order to avoid future dependencies:
  - libbluetooth3 as library to use the bluez linux bluetooth stack.

Our libraries has been developed with **4.96 version** (bluez 4.96-0ubuntu4 package) of BlueZ (latest release available for download is  $4.101^7$  nowadays). Upper versions of 4.80 is mandatory due to it was the first release where HDP was integrated.

API of GNU/Linux official bluetooth stack, BlueZ, uses the popular

<sup>7 &</sup>lt;a href="http://www.bluez.org/download/">http://www.bluez.org/download/</a>

message bus system created by Red Hat company D-Bus<sup>8</sup>. D-Bus is an inter-process communication (IPC) open-source system that enables a reliable way for exchanging data (messages) directly between software applications and between software applications and the operating system (where the operating system would typically include the kernel and any system daemons or processes) over the same desktop session. Although there are some implementations<sup>9</sup> (bindings) of D-Bus in several higher level development languages (Java, Python, Qt, Glib, C#, etc) we decided to base our development over the low level D-Bus implementation which API is written in C code for the following reasons:

- Low level implementation has no additional dependencies unlike higher level bindings and it is distributed by default in GNU/Linux OS also.
- Java binding (the closest to uAAL partners interests, for obvious reasons) has a very important lack compared with low level implementation (C code): file descriptor data type is not supported in the latest release available for download (D-Bus works with Unix sockets so a valid file descriptor should be reserved to any HDP channel created prior to send and receive data bytes between agents and managers).

Our libraries has been developed with **1.4.14 D-Bus version** (latest release available for download is 1.6.1 nowadays). Please keep in mind that working with low level D-Bus C API is not easy and has a very pronounced learning curve so you should only go deeper just if necessary (main page of D-Bus C API<sup>10</sup> welcomes visitors with a warm: *If you use this low-level API directly, you're signing up for some pain*). Please, ensure that any version of D-Bus is conveniently installed. Additionally, as developers we need to install some packages related with **libdbus** library in order to avoid future dependencies:

- libglib2.0-dev as development files for GLib libraries.
- libdbus1.0-cil-dev as CLI implementation of D-Bus.
- libdbus-glib1.0-cil-dev as CLI implementation of D-Bus (Glib mainloop integration).
- libdbus-1-dev as development headers for simple interprocess messaging system.
- libdbus-glib-1-dev as Glib interface.

<sup>8 &</sup>lt;a href="http://www.freedesktop.org/wiki/Software/dbus">http://www.freedesktop.org/wiki/Software/dbus</a>

<sup>9 &</sup>lt;a href="http://www.freedesktop.org/wiki/Software/DBusBindings">http://www.freedesktop.org/wiki/Software/DBusBindings</a>

<sup>10</sup> http://dbus.freedesktop.org/doc/api/html/index.html

Е	Paquete	Versión instalada	Última versión	Descripción
		0.5.0-3build1	0.5.0-3build1	CLI implementation of D-Bus (GLib mainloop integration)
		0.5.0-3build1	0.5.0-3build1	CLI implementation of D-Bus (GLib mainloop integration) - development files
		0.7.0-4	0.7.0-4	CLI implementation of D-Bus
	libdbus1.0-cil-dev	0.7.0-4	0.7.0-4	Implementación CLI de D-Bus - archivos de desarrollo
		1.4.14-1ubuntu1	1.4.14-1ubuntu1	simple interprocess messaging system (development headers)
		0.94-4	0.94-4	sistema de mensajes interproceso simple (interfaz GLib)
		0.94-4	0.94-4	sistema simple de mensajería inter-proceso (biblioteca dinámica basada en GLib)
	dbus      d	1.4.14-1ubuntu1	1.4.14-1ubuntu1	simple interprocess messaging system (daemon and utilities)
	ibdbus-1-3	1.4.14-1ubuntu1	1.4.14-1ubuntu1	simple interprocess messaging system (library)

Figure 4: Synaptic package manager view with "libdbus" and "installed packages" as quick filters

If everything goes as expected, our local folders (usr/lib/i386-linux-gnu/dbus-1.0/include/dbus and usr/include/dbus-1.0/dbus) Should contain some new header files as shown following:

```
🚫 🖨 📵 angel@angel-laptop: ~
angel@angel-laptop:~$ ls -l /usr/lib/i386-linux-gnu/dbus-1.0/include/dbus/
total 4
-rw-r--r-- 1 root root 2089 2011-09-02 11:30 dbus-arch-deps.h
angel@angel-laptop:~$ ls -l /usr/include/dbus-1.0/dbus/
total 204
-rw-r--r-- 1 root root 2135 2011-09-02 11:30 dbus-address.h
-rw-r--r-- 1 root root 3472 2011-09-02 11:30 dbus-bus.h
-rw-r--r-- 1 root root 25953 2012-04-10 17:02 dbus-connection.h
rw-r--r-- 1 root root 2886 2011-09-02 11:30 dbus-errors.h
 rw-r--r-- 1 root root 23000 2011-07-07 16:52 dbus-glib-bindings.h
 rw-r--r-- 1 root root 14376 2011-07-07 16:51 dbus-glib.h
 rw-r--r-- 1 root root 2875 2011-07-07 16:51 dbus-glib-lowlevel.h
rw-r--r-- 1 root root 1384 2011-07-07 16:51 dbus-gvalue-parse-variant.h
-rw-r--r-- 1 root root 3925 2012-04-10 17:03 dbus.h
-rw-r--r-- 1 root root 4881 2011-09-02 11:30 dbus-macros.h
-rw-r--r-- 1 root root 1975 2011-09-02 11:30 dbus-memory.h
-rw-r--r-- 1 root root 13049 2011-09-02 11:30 dbus-message.h
-rw-r--r-- 1 root root 1609 2011-09-02 11:30 dbus-misc.h
-rw-r--r-- 1 root root 23178 2011-09-02 11:30 dbus-protocol.h
-rw-r--r-- 1 root root 3534 2011-07-27 13:30 dbus-python.h
rw-r--r-- 1 root root 4794 2011-09-02 11:30 dbus-server.h
 rw-r--r-- 1 root root 5135 2011-09-02 11:30 dbus-shared.h
-rw-r--r-- 1 root root 2956 2011-09-02 11:30 dbus-signature.h
-rw-r--r-- 1 root root 9004 2011-09-02 11:30 dbus-threads.h
-rw-r--r-- 1 root root 3635 2012-04-10 17:10 dbus-types.h
angel@angel-laptop:~$
```

Figure 5: D-Bus IPC system header files

At this point, we should ensure that D-Bus, BlueZ and our bluetooth dongle works fine as a whole. In order to verify (again) that your bluetooth adapter is ready to work run hciconfig command (for example) as root.

```
angel@angel-laptop:~

angel@angel-laptop:~$ hciconfig
hci0: Type: BR/EDR Bus: USB
BD Address: 00:16:38:C1:9D:EA ACL MTU: 1017:8 SCO MTU: 64:0
UP RUNNING PSCAN ISCAN
RX bytes:1005 acl:0 sco:0 events:27 errors:0
TX bytes:622 acl:0 sco:0 commands:27 errors:0

angel@angel-laptop:~$
```

Figure 6: Checking availability of local bluetooth adapter

Next, type the following D-Bus command to get your assigned BlueZ D-Bus address:

```
$ dbus-send --system --type=method_call --print-reply
--dest=org.bluez "/" org.bluez.Manager.ListAdapters
```

**Table 3: List of local bluetooth adapters** 

```
    angel@angel-laptop: ~
angel@angel-laptop: ~
before the discrete angel@angel-laptop: ~
angel@angel-laptop: ~
before the discrete angel@angel-lapto
```

Figure 7: Local bluetooth adapter path returned by D-Bus system

The returned object path as string org/bluez/1043/hci0 is our adapter's name where 1043 is the PID (Process ID) bluetooth process at this session (Important note: local bluetooth adapter PID varies after each computer boot so it shouldn't be consider as static value). As next step, we can (for instance) request for local bluetooth adapter properties as follows:

```
$ dbus-send --system --type=method_call --print-reply
--dest=org.bluez "/org/bluez/1043/hci0"
org.bluez.Adapter.GetProperties
```

Table 4: Default local bluetooth adapter properties

If everything goes well, automatically output should be (more or less) similar to:

```
method return sender=:1.8 -> dest=:1.74 reply_serial=2
  array [
     dict entry(
        string "Address"
                           string "00:16:38:C1:9D:EA"
        variant
     dict entry(
        string "Name"
                           string "angel-laptop-0"
        variant
     dict entry(
        string "Class"
        variant
                          uint32 4849920
     dict entry(
        string "Powered"
                          boolean true
        variant
     dict entry(
        string "Discoverable"
        variant
                   boolean true
     dict entry(
       string "Pairable"
        variant
                          boolean true
     dict entry(
        string "DiscoverableTimeout"
        variant
     dict entry(
        string "PairableTimeout"
        variant
                          uint32 0
     dict entry(
        string "Discovering"
        variant boolean false
     dict entry(
        string "Devices"
                            array [
              object path "/org/bluez/1043/hci0/dev_00_09_1F_80_0A_E0"
              object path "/org/bluez/1043/hci0/dev 00 09 1F 80 04 D6"
     dict entry(
        string "UUIDs"
        variant
                            array [
              string "00001000-0000-1000-8000-00805f9b34fb"
              string "00001001-0000-1000-8000-00805f9b34fb"
              string "00001112-0000-1000-8000-00805f9b34fb"
              string "0000111f-0000-1000-8000-00805f9b34fb"
              string "0000110a-0000-1000-8000-00805f9b34fb"
              string "0000110c-0000-1000-8000-00805f9b34fb"
              string "0000110e-0000-1000-8000-00805f9b34fb"
              string "00001103-0000-1000-8000-00805f9b34fb"
```

Table 5: Default local bluetooth adapter properties returned by D-Bus system

#### Shoulders to the wheel

Summarizing, for being able to run the libraries, you will have to need a PC with GNU/Linux operating system installed with one of the latest versions of BlueZ library and D-Bus system for managing Bluetooth connections. In our case, we had a 32-bit Ubuntu 11.10 with BlueZ 4.96 and Oracle Java JDK 7 also.

<u>Important note</u>: as prior step to exchange data bytes between Continua agents and managers both should have been paired first. Developed libraries doesn't have an automatic process to manage this task at this release so you should use any bluetooth manager software available for download. We decided to use the software distributed with the OS by default for simplicity. Next figure shows that application with both the weigh scale and the upper arm blood pressure monitor paired with our computer (you will probably need to consult the right instruction manual of chosen Continua device to get the secret PIN required as final step of any pairing process).




Figure 8: Bluetooth software manager with Continua compliant paired devices

Once the Continua agents and our PC (or any smart device with manager role) has been paired, we can continuous with the next steps as detailed below.

#### **Working without uAAL middleware**

As uAAL middleware runs inside an OSGi framework and this environment is not too intuitive and easy to integrate third party libraries (at the beginning), a good practice could be start working with the libraries developed in standalone mode as first step. Please, follow next steps:

- 1. Update your local copy of the uAAL repository with the latest version of the libraries. The contents you must have should be:
  - Jar file with the implementation of ISO/IEEE 11073 standard (java pure code, distributed in several packages and without external dependencies): x73implementation.jar
  - Native library written in C code (with ".so" extension) with JNI functions for accessing Bluetooth stack with Java through D-Bus IPC system: libHDPnative NOUAAL.so
  - Javadoc documentation for helping uAAL developers creating Continua managers with Bluetooth HDP capabilities: javadoc.rar
  - Java code (Eclipse IDE project format) with a manager as sample: BloodPressureMonitor\_x073\_Sample

With this three elements and the sample of code, you will be able to create your own managers, and obtain measurements from any Continua compliant medical device such as a weighing scale or a blood pressure monitor.

- 2. Ensure that your Continua device/devices is/are successfully paired with your computer (this step is mandatory, don't forget it).
- 3. Create a new folder **exactly named** jni\_lib at your home path (Important note: this path must be exact. Otherwise, it won't work).
- 4. Copy the libHDPnative\_NOUAAL.so file into the previously created directory.

```
angel@angel-laptop:~/jni_lib

angel@angel-laptop:~/jni_lib$ ls -l

total 64

-rwxrwxr-x 1 angel angel 31197 2012-06-29 10:07 libHDPnative_NOUAAL.so

-rwxrwxr-x 1 angel angel 32044 2012-06-27 13:15 libHDPnative.so

angel@angel-laptop:~/jni_lib$
```

Figure 9: Native library inside our local folder

5. Import the sample project (BloodPressureMonitor\_x073\_Sample) available in uAAL repository in Eclipse IDE or create a new one with the same package structure than the other one. In any case, jar package with ISO/IEEE x073 implementation must be imported into the Java project. So, your software project should contain next classes/imported packages at least:

- hdpManager.java is the core of the manager application. It should implement hdpManagerListener.java interface.
- hdpManagerListener.java is an interface that only defines three selfexplanatory methods: onChannelConnected, onChannelDisconnected and onDataReceived. Please, refer the aforementioned Javadoc file to get more detailed information.
- x73implementation.jar is the ISO/IEEE 11073 standard implementation in Java code. It is responsible for processing and managing x073 received data frames.




Figure 10: Package explorer of Eclipse IDE

Important note: package name that contains hdpManager.java and hdpManagerListener.java classes must be bluetooth. Otherwise, native library will return an error (this file is, more or less, autogenerated by the javah tool included in the JDK so this value is the right expected).

- 6. Since the output text messages from the C functions at libHDPnative.so file seems not to be showing in the console of the Eclipse IDE platform, we recommend to execute the software from the GNU/Linux OS console terminal as follows:
  - In hdpManager.java file change the macAddressRemoteDevice attribute with the corresponding MAC value (line 76) of your blood pressure device.
  - If your project is built automatically by Eclipse IDE (Project → Build automatically option checked) go to the bin folder inside the sample project directory (workspace) and execute java

bluetooth.hdpManager command (where hdpManager is the manager main class and bluetooth the name of the container package). Next table shows some properties about the local bluetooth adapter, its valid PID and how an HDP channel is waiting for Continua agent frames to start the exchange data process.

```
angel@angel-laptop: ~/TSB/src/workspace_uaal/BloodPressureMonitor_x073_Sample/bin
angel@angel-laptop:~/TSB/src/workspace_uaal/BloodPressureMonitor_x073_Sample/bin$ java bluetooth.hdpManager
UNASSOCIATED
Default bluetooth adapter path: /org/bluez/1107/hci0
List of local bluetooth adapters:
Adapters: /org/bluez/1107/hci0
Local bluetooth properties:
Address: 00:16:38:C1:9D:EA
Name: angel-laptop-0
Class: 4849920
Powered: TRUE
Discoverable: TRUE
Pairable: TRUE
DiscoverableTimeout: 0
PairableTimeout: 0
Discovering: FALSE
Devices: /org/bluez/1107/hci0/dev_00_09_1F_80_0A_E0
/org/bluez/1107/hci0/dev_00_09_1F_80_04_D6
UUIDs: 00001000-0000-1000-8000-00805f9b34fb
00001001-0000-1000-8000-00805f9b34fb
00001112-0000-1000-8000-00805f9b34fb
0000111f-0000-1000-8000-00805f9b34fb
0000110a-0000-1000-8000-00805f9b34fb
0000110c-0000-1000-8000-00805f9b34fb
0000110e-0000-1000-8000-00805f9b34fb
```

Figure 11: Blood pressure manager successfully launched

7. Start getting measurements from the Continua devices in our area of coverage is the actual step. Next images shows an example of a blood pressure measurement:

```
LOG: Initiating processing of a received measurement from a Blood Pressure
MEASUREMENT: 114 mmHg
MEASUREMENT: 65 mmHg
MEASUREMENT: 83 mmHg
MEASUREMENT: 87 bpm
MEASUREMENT: Date of measurement: 11:49:05 - 17/07/2012
SYSTEM TIME: 11:49:07 - 17/07/2012
```



Figure 12: Measurements done with a real Continua blood pressure monitor

```
$ pwd
/.../workspace uaal/BloodPressureMonitor x073 Sample/bin
/.../workspace uaal/BloodPressureMonitor x073 Sample/bin$ java bluetooth.hdpManager
UNASSOCIATED
Default bluetooth adapter path: /org/bluez/1107/hci0
List of local bluetooth adapters:
Adapters: /org/bluez/1107/hci0
Local bluetooth properties:
Address: 00:16:38:C1:9D:EA
Name: angel-laptop-0
Class: 4849920
Powered: TRUE
Discoverable: TRUE
Pairable: TRUE
DiscoverableTimeout: 0
PairableTimeout: 0
Discovering: FALSE
Devices: /org/bluez/1107/hci0/dev 00 09 1F 80 0A E0
/org/bluez/1107/hci0/dev 00 09 1F 80 04 D6
UUIDs: 00001000-0000-1000-8000-00805f9b34fb
00001001-0000-1000-8000-00805f9b34fb
00001112-0000-1000-8000-00805f9b34fb
0000111f-0000-1000-8000-00805f9b34fb
0000110a-0000-1000-8000-00805f9b34fb
0000110c-0000-1000-8000-00805f9b34fb
0000110e-0000-1000-8000-00805f9b34fb
00001103-0000-1000-8000-00805f9b34fb
Remote bluetooth adapter path: /org/bluez/1107/hci0/dev_00_09_1F_80_04_D6
Remote device properties (BEFORE):
Address: 00:09:1F:80:04:D6
Name: AND BP5091100122
Alias: AND BP5091100122
Class: 2308
Paired: TRUE
Trusted: TRUE
Blocked: FALSE
Connected: FALSE
UUIDs: 00001200-0000-1000-8000-00805f9b34fb
00001401-0000-1000-8000-00805f9b34fb
Services:
Adapter: /org/bluez/1107/hci0
Property Trusted successfully changed
Remote device properties (AFTER):
Address: 00:09:1F:80:04:D6
Name: AND BP5091100122
Alias: AND BP5091100122
Class: 2308
Paired: TRUE
Trusted: TRUE
Blocked: FALSE
Connected: FALSE
UUIDs: 00001200-0000-1000-8000-00805f9b34fb
00001401-0000-1000-8000-00805f9b34fb
Services:
Adapter: /org/bluez/1107/hci0
HDP application identifier: /org/bluez/health app 3
HDP application (/org/bluez/health app 3) successfully closed
HDP application identifier: /org/bluez/health_app_4
DBUS sytem status: Up
CONNECTED UNASSOCIATED
Waiting for source connections...
```

Table 6: Continua manager waiting for input source connections (agents)

#### **Working with uAAL middleware**

Once we were able to get real time measurements from Continua devices in standalone mode as has shown before we need to migrate that technological solution inside OSGi environments due to this framework is the basis of uAAL middleware. Please, follow next steps and be patient:

- 1. Update your local copy of the uAAL repository with the latest version of the libraries. The contents you must have should be:
  - Maven project with the implementation of ISO/IEEE 11073 standard (java pure code, distributed in several packages and without external dependencies): ieee.x73.std
  - Native library written in C code (with ".so" extension) with JNI functions for accessing Bluetooth stack with Java through D-Bus IPC system (please, choose the most appropriate version according your OS): libHDPnative.so
  - Javadoc documentation for helping uAAL developers creating Continua managers with Bluetooth HDP capabilities: javadoc.rar
  - Some Maven projects as code samples:
    - 1. weighingScalePublisher x073 Sample: Simple app without GUI.
    - 2. ContinuaManager: advanced app with GUI enabled.
    - 3. "UniversAALized" sample: more advanced app following uAAL design patterns. It is composed by four components:
      - continua manager ontology: ont.continuahealthmanager
      - continua manager service: continua.hdp.manager
      - continua manager caller: continua.hdp.manager.caller.test
      - continua manager subscriber:
         continua.hdp.manager.subscriber.test

With this three elements and the samples of code, you will be able to create your own managers, obtain measurements from any Continua compliant medical device such as a weighing scales or a blood pressure monitors and send data received to uAAL context bus.

- 2. Ensure that your Continua device/devices is/are successfully paired with your computer (this step is mandatory, don't forget it).
- 3. Create a new folder **exactly named** jni\_lib at your home path (Important note: this path must be exact. Otherwise, it won't work).
- 4. Copy the libHDPnative.so file into the previously created directory (please, review figure 9).
- 5. Import one of the sample Maven projects from your local copy of the uAAL repository to your Eclipse IDE. In order to work with our more realistic sample we need to import four Maven projects:

 File → Import → Existing Maven Projects → continua.hdp.manager Import Existing Projects Maven continua.hdp.manager.caller.test File Import Maven Projects Existing continua.hdp.manager.subscriber.test} ▽ □ 🛮 🗏 Package Explorer 🔀 continua.hdp.manager [sandboxes/amartinez/continua.hdp.manager] \* continua.hdp.manager.caller.test [sandboxes/amartinez/continua.hdp.manager.caller.test] continua.hdp.manager.subscriber.test [sandboxes/amartinez/continua.hdp.manager.subscriber.test] ContinuaManager [ContinuaManager/ContinuaManager] ContinuaManagerClientGui v0 ContinuaManagerClientGui v1 continuaWeightingScalePublisher\_v0 continuaWeightScalePublisher i drools.publisher.test drools.suscriber.test [sandboxes/mllorente/drools.suscriber.test] libraries/ieee.x73.std | icee.x73.std | icee.x73.st ieee.x73.std\_backup\_29\_06\_2012 iddi.hw.exporter.personalhealthdevice.simulator [sandboxes/thomas.fuxreiter/lddi.hw.exporter.per Itba.manager [sandboxes/mllorente/ltba.manager] Similar ont.continuahealthmanager [sandboxes/amartinez/ont.continuahealthmanager] org.universAAL.ontology.x73 i samples.pom mp.lighting.client.osgi mp.lighting.server.osgi mp.lighting.uiclient 

File → Import → Existing Maven Projects → ont.continuahealthmanager

What should be expected by each Maven project imported:

- Continua.hdp.manager is a Continua Health manager service working as a whole (you could switch on and off it). It is the main responsible for managing secure connections between agents and managers and for processing low level frames. A new independent and isolated Continua manager is created when an uAAL caller message is received.
- Continua.hdp.manager.caller.test is an uAAL caller class. Just for testing.
- Continua.hdp.manager.subscriber.test is an uAAL subscriber class. Just for testing.
- Ont.continuahealthmanager is the Continua Health Alliance Manager

ontology. It is used by the Continua Health manager service.

Important note: (1) don't change packages names. Otherwise, native library will return an error (this file is, more or less, auto-generated by the javah tool included in the JDK so this value is the right expected. In conclusion: if thoses package are renamed, functions in native library should be renamed too. Otherwise, both names will not match and an error will be return as output) (2) manage possible external dependencies with the pom file (pom.xml). There are some ontologies specifically created to model Continua agents and managers. Please, check them before launching the service. as weighing scales and blood pressures. Please review org.univershal.ontology.x73 package in uAAL SVN repositories (3) Delete (if proceed) the x073implementation.jar JAR file previously copied inside your local path folder (/usr/lib/jvm/java-7-oracle/jre/lib/ext) to force maven project get the library from the right source.

6. In GUI.java file (Continua.hdp.manager.caller.test project) you should change the remoteMacAddress attribute with the corresponding MAC value (lines 189 and 216) of your weighing scale device and blood pressure monitor.

```
J GUI.java 器
 180 ** Action listener for radio buttons and jbuttons */
△181⊖ublic void actionPerformed(ActionEvent e) {
        if(e.getActionCommand().equals("bloodPressure")) {
 182
 183
            // Blood pressure device
 184
            if(realMeasurementButton.isSelected()) {
 185⊖
                new Thread(){
▲186⊝
                    public void run() {
                         realMeasurement = true;
                         remoteMacAddress = "00:09:1F:80:04:D6";
 189
                         simulatedMeasurementButton.setEnabled(false);
 191
                         System.out.println("Continua Health Manager service request...");
 192
 193
                         caller = new DefaultServiceCaller(mcx);
 194
                         ServiceRequest addReq = new ServiceRequest(new ContinuaHealthManager(null), null);
                         addReq.addAddEffect(new String[] { ContinuaHealthManager.PROP_HAS_CONTINUA_DEVICE },
 195
2196
                                 //TODO Put here the right MAC address of your Continua device
                                 new ContinuaHealthDevice(ContinuaHealthManagerOntology.NAMESPACE + "new Continua
 197
 198
                         System.out.println(addReq.toString());
```

6. Launch this application with an appropriate Run Configurations. Please, check that all previously installed bundles are Active less than our Caller bundle. (remember: output text messages from the C functions at libHDPnative.so file seems not to be showing in the console of the Eclipse IDE platform, this bug will be solved ASAP by TSB team):

```
-> DS
START LEVEL 20
      State
                     Level Name
  ID
   0] [Active
                   ] [
                         0] System Bundle (2.0.1)
   1] [Active
                         8] universAAL ontology for personal health devices (1.2.1.SNAPSHOT)
   2] [Active
                   ] [
                         8] uAAL Continua Manager service (1.2.2.SNAPSHOT)
   3] [Active
                   ] [
                         8] uAAL Continua Manager service subscriber test (1.2.2.SNAPSHOT)
                         /] The Physical World Unitology (1.2.1.SNAPSHUT)
   4] [ACTIVE
                         7] uAAL ontology for Continua Health Alliance managers (1.2.2.SNAPSHOT)
   5] [Active
                   ] [
   oj [ACTIVE
                         of the measurement untology (1.2.0)
   7] [Installed
                         9] uAAL Continua Manager service caller test (1.2.2.SNAPSHOT)
   o] [ACTIVE
                         oj user interaction bus (usul) (1.2.1.5NAPSHUT)
   9] [Active
                         5] Context Bus (OSGi) (1.2.1.SNAPSHOT)
                   ] [
  10] [Active
                         5] Service Bus (OSGi) (1.2.1.SNAPSHOT)
                   ] [
  11] [Active
                         5] Unit (1.2.0)
  12] [Active
                         5] wrap mvn org.universAAL.continua ieee.x73.std 1.1.0 (0)
  13] [Active
                   1 [
                         4] Data Representation (OSGi) (1.2.1.SNAPSHOT)
  14] [Active
                         4] RDF/OWL Turtle serializer (OSGi) (1.2.1.SNAPSHOT)
                   ] [
                         3] Bus Model (OSGi) (1.2.1.SNAPSHOT)
  15] [Active
                   ] [
  16] [Active
                         2] OSGi R4 Core Bundle (4.1)
                   ] [
  17] [Active
                         2] Apache Felix Configuration Admin Service (1.2.4)
                   ] [
  18] [Active
                   ] [
                         2] Container Interfaces (OSGi) (1.2.1.SNAPSHOT)
  19] [Active
                         2] wrap mvn jp.go.ipa jgcl 1.0 (0)
  20] [Active
                         2] Sun Microsystems' JDBC Driver for MySQL (5.1.17)
                   1 [
                         2] OPS4J Pax Logging - API (1.6.2)
  21] [Active
                   ] [
  22] [Active
                   ] [
                         2] wrap mvn java3d vecmath 1.3.1 (0)
  23] [Active
                   ] [
                         2] ACL Interfaces (OSGi) (1.2.1.SNAPSHOT)
  24] [Active
                         2] wrap_mvn_org.bouncycastle_jce.jdk13_144 (0)
                   ] [
  25] [Active
                   ] [
                         2] wrap_mvn_java3d_j3d-core_1.3.1 (0)
  26] [Active
                         2] wrap mvn org.osgi osgi R4 compendium 1.0 (0)
  27] [Active
                         2] OSGi Container (1.2.1.SNAPSHOT)
  28] [Active
                   ] [
                         2] Jena OSGi Bundle (1.0.0)
                         1] Apache Felix Shell Service (1.4.1)
  29] [Active
                   ] [
[ 30] [Active

    Apache Felix Shell TUI (1.4.1)
```

<u>Important note</u>: <u>VM arguments</u> at <u>Arguments</u> tab in the <u>Run</u><u>Configurations</u> window should contain the path where native library
<u>has been stored as follows: -Djava.library.path=/home/user/jni lib</u>

```
-Dosgi.noShutdown=true -Dfelix.log.level=4
-Dorg.universAAL.middleware.peer.is_coordinator=true-
-Dorg.universAAL.middleware.peer.member_of=urn:org.universAAL.aal_space:test_env-
-Dbundles.configuration.location=$(workspace_loc)/rundir/confadmin-
-Djava.library.path=/home/angel/jni_lib
```

7. Start the Caller bundle with the appropriate ID number. If all goes as expected a simple GUI should be shown in your screen.



- 8. If you have access to real Continua devices, press the appropriate combination of checkbox/button and wait for agents events. If not, you can send simulated and randomized measures to the uAAL context bus. At any case, when the measurement is published in the context bus, the subscriber is able to get the data and show it in another application.
  - 1. Subscriber for simulated measurements (blood pressure monitor, for instance).



```
@prefix rdf: <http://www.w3.org/1999/02/22-rdf-syntax-ns#> .
@prefix ns3: <http://ontology.universAAL.org/Device.owl#> .
@prefix : <http://http://ontology.universAAL.org/PersonalHealtDeviceSimulator.owl#>
<urn:org.universAAL.middleware.context.rdf:ContextEvent#13b8f7952e8@angel-laptop+13fbf</pre>
 a ns:ContextEvent;
  rdf:subject :BloodPressureMonitor ;
 ns:hasTimestamp "1355321805489"^^xsd:long;
  rdf:predicate ns2:hasMeasurement;
  rdf:object :BloodPressureMeasurement .
ns:gauge a ns:ContextProviderType .
:BloodPressureMeasurement ns2:hasMeasuredBPSys :Measurement Systolic ;
 a ns2:BloodPressureMeasurement ,
    ns3:Device ,
    pvn:PhysicalThing;
 ns2:hasMeasuredBPDia :Measurement Diastolic ;
 ns2:hasMeasuredHeartRate :Measurement Heartrate
:Measurement Systolic a ns1:Measurement ;
 ns1:value "114"^^xsd:string .
:PersonalHealthDeviceContextProvider a ns:ContextProvider ;
 ns:hasType ns:gauge .
:Measurement Heartrate a ns1:Measurement ;
 ns1:value "149"^^xsd:string .
:Measurement Diastolic a ns1:Measurement ;
 ns1:value "71"^^xsd:string .
:BloodPressureMonitor ns2:hasMeasurement :BloodPressureMeasurement ;
 a ns2:BloodPressureMonitor ,
   ns3:Device ,
   pvn:PhysicalThing .
```

2. Subscriber for real measurements (weighing scale, for instance).



```
@prefix ns1: <http://ontology.universaal.org/Measurement.owl#> .
@prefix xsd: <http://www.w3.org/2001/XMLSchema#> .
@prefix ns2: <http://ontology.universAAL.org/PersonalHealthDevice.owl#> .
@prefix rdf: <http://www.w3.org/1999/02/22-rdf-syntax-ns#> .
@prefix : <http://ontology.universAAL.org/Context.owl#> .
<urn:org.universAAL.middleware.context.rdf:ContextEvent#13b8f7952e8@angel</pre>
  a :ContextEvent :
  rdf:subject ns:WeighingScale;
  :hasTimestamp "1355322050968"^^xsd:long ;
  rdf:predicate ns2:hasMeasuredWeight;
  rdf:object ns:Measurement Weight .
:gauge a :ContextProviderType .
ns:PersonalHealthDeviceContextProvider a :ContextProvider ;
  :hasType :gauge .
ns:WeighingScale a ns2:WeighingScale ,
    <http://ontology.universAAL.org/Device.owl#Device> ,
    <http://ontology.universAAL.org/uAAL.owl#PhysicalThing> ;
  ns2:hasMeasuredWeight ns:Measurement Weight .
ns:Measurement Weight a ns1:Measurement ;
  ns1:value "3500"^^xsd:string .
```

<u>Important note</u>: When the Caller application is closed an <u>Invalid</u> <u>BundleContext</u> <u>exception</u> is thrown. It seems that is a well known problem inside uAAL developers community so we assume it as normal response. Any advice, suggestion to avoid it would be appreciated.

```
Exception in thread "Thread-14" java.lang.IllegalStateException: Invalid BundleContext.

at org.apache.felix.framework.BundleContextImpl.checkValidity(BundleContextImpl.java:393)

at org.apache.felix.framework.BundleContextImpl.getBundle(BundleContextImpl.java:92)

at org.universAAL.middleware.container.osgi.uAALBundleContext.getID(uAALBundleContext.java:118)

at org.universAAL.middleware.container.utils.LogUtils.log(LogUtils.java:105)

at org.universAAL.middleware.service.ServiceCaller.handleReply(ServiceCaller.java:163)

at org.universAAL.middleware.service.impl.ServiceStrategy.replyToLocalCaller(ServiceStrategy.java:1584)

at org.universAAL.middleware.service.impl.ServiceStrategy.sendServiceResponse(ServiceStrategy.java:761)

at org.universAAL.middleware.service.impl.ServiceStrategy.handle(ServiceStrategy.java:1013)

at org.universAAL.middleware.sodapop.BusStrategy$HandlerThread.run(BusStrategy.java:50)
```

#### **Troubleshooting**

1. <u>Error</u>: compiling any source class always returns a java.lang.UnsupportedClassVersion. <u>Possible solution</u>: Change your JDK to the official Oracle Java JDK or update its version.

```
$ java -version
java version "1.7.0_04"
Java(TM) SE Runtime Environment (build 1.7.0_04-b20)
Java HotSpot(TM) Server VM (build 23.0-b21, mixed mode)
```

Table 7: UnsupportedClassVersion error

2. <u>Error</u>: launching manager application with x073implementation.jar file imported always returns a java.lang.NoClassDefFoundError. <u>Possible solution</u>: Check the classpath or directly copy the jar file inside /usr/lib/jvm/java-7-oracle/jre/lib/ext path folder.

```
🙆 同 🏻 angel@angel-laptop: ~/TSB/src/workspace_uaal/BloodPressureMonitor_x073_Sample/bin
angel@angel-laptop:~/TSB/src/workspace_uaal/BloodPressureMonitor_x073_Sample/bin$ java bluetooth.hdpManager
Exception in thread "main" java.lang.NoClassDefFoundError: events/EventIEEEListener
           at java.lang.ClassLoader.defineClass1(Native Method)
at java.lang.ClassLoader.defineClass(ClassLoader.java:791)
           at java.security.SecureClassLoader.defineClass(SecureClassLoader.java:142) at java.net.URLClassLoader.defineClass(URLClassLoader.java:449)
           at java.net.URLClassLoader.access$100(URLClassLoader.java:71)
           at java.net.URLClassLoader$1.run(URLClassLoader.java:361)
           at java.net.URLClassLoader$1.run(URLClassLoader.java:355)
           at java.security.AccessController.doPrivileged(Native Method)
           at java.net.URLĆlassLoader.findClass(URLClassLoader.java:354)
at java.lang.ClassLoader.loadClass(ClassLoader.java:423)
           at sun.misc.LauncherSAppClassLoader.loadClass(Launcher.java:308) at java.lang.ClassLoader.loadClass(ClassLoader.java:356) at sun.launcher.LauncherHelper.checkAndLoadMain(LauncherHelper.java:480)
Caused by: java.lang.ClassNotFoundException: events.EventIEEEListener at java.net.URLClassLoader$1.run(URLClassLoader.java:366) at java.net.URLClassLoader$1.run(URLClassLoader.java:355) at java.security.AccessController.dClassLoader(Native Method)
           at java.net.URLClassLoader.findClass(URLClassLoader.java:354)
           at java.lang.ClassLoader.loadClass(ClassLoader.java:423)
           at sun.misc.Launcher$AppClassLoader.loadClass(Launcher.java:308)
           at java.lang.ClassLoader.loadClass(ClassLoader.java:356)
                 13 more
```

Figure 13: NoClassDefFoundError error

3. <u>Error</u>: launching manager application always returns an undefined symbol error. <u>Possible solution</u>: Execute 1dd -d command (print shared library dependencies) in the OS console in the folder where native libraries (remember: .so files) where stored. If these libraries has successfully been linked in your system, you should get an output similar as follows:

Figure 14: Native libraries successfully linked

4. <u>Error</u>: launching manager application always returns a java.lang.UnsatisfiedLinkError. <u>Possible solution</u>: Check the classpath or update your LD LIBRARY PATH global variable.

#### **Appendix A**

#### Setting up uAAL development environment on GNU/Linux

Important note: this guidelines are not official documentation supported and maintained by universAAL Consortium. This document has been elaborated (of our own necessity) in order to help any GNU/Linux user starting setting up the minimum required environment to run and develop AAL services based on uAAL middleware. Official information available in the universAAL Developer Handbook<sup>11</sup> only covers this process in Microsoft Windows environments so due to this lack we decide to create this simple guide from scratch. As any first release it isn't exempt from errors, so we apologize for any inconvenience. Please, feel free to comment or modify anything of this appendix. Any suggestion will be appreciated.

We will assume that you are using a 32-bit Linux Ubuntu 11.10 operating system.

- 1. Install an appropriate JDK. We recommend Oracle-jdk7.
- 2. Create directory tree. For example, we will use as root of universaal environment a folder in home directory.
  - 1. Create ~/universAAL folder. Please note that ~ symbol refers to home directory.
- 3. Install maven software:
  - 1. Download the latest version package from the official URL<sup>12</sup>
  - 2. Extract the files in ~/universAAL/maven. This will as seen in screenshot:



3. Set up the proper environment variables. For doing this permanent, we will have to edit the .bashrc file in home folder. In console terminal,

<sup>11</sup> http://forge.universaal.org/gf/

<sup>12</sup> http://maven.apache.org/download.html

#### write gedit ~/.bashrc (at the end of that text file, add next lines):

```
export M2_HOME=/home/username/universAAL/maven/apache-maven-3.0.4
export M2=$M2_HOME/bin
export MAVEN_OPTS="-Xms256m -Xmx512m"
export PATH=$M2:$PATH
export JAVA_HOME=/usr/lib/jvm/java-7-oracle
export PATH=$JAVA_HOME/bin:$PATH
```

4. Check if it is working correctly. Enter the mvn --version in console terminal. If everything is okay, we have to see information about our computer and the version of maven software just installed:

```
angel@angel-laptop:~$ mvn --version
Apache Maven 3.0.4 (r1232337; 2012-01-17 09:44:56+0100)
Maven home: /home/angel/TSB/uAAL/Maven/apache-maven-3.0.4
Java version: 1.7.0_04, vendor: Oracle Corporation
Java home: /usr/lib/jvm/java-7-oracle/jre
Default locale: es_ES, platform encoding: UTF-8
OS name: "linux", version: "3.0.0-19-generic-pae", arch: "i386", family: "unix"
angel@angel-laptop:~$
```

- 5. Once we have maven software running in our machine, it is time to install a SVN client and to obtain the universAAL resources (source code). This guidelines will assume that you are registered in universAAL GForce web page<sup>13</sup> and you have been granted access to Middleware, Ontologies and Support projects.
  - Install any SVN client which can be retrieved from Ubuntu Software Center or similar (for instance: RabbitVCS<sup>14</sup> or SVN Workbech<sup>15</sup>).
  - Create the folders (for instance):
    - o ~/universAAL/svn/middleware/
    - o ~/universAAL/svn/ontologies/
    - o ~/universAAL/svn/support/
  - Get uAAL source code from universAAL repository accordingly with SVN client previously chosen. Wait until all files are downloaded.
    - SVN source: http://forge.universaal.org/svn/middleware
    - SVN source: http://forge.universaal.org/svn/ontologies
    - SVN source: http://forge.universaal.org/svn/support

<sup>13 &</sup>lt;a href="http://forge.universaal.org/gf/">http://forge.universaal.org/gf/</a>

<sup>14 &</sup>lt;a href="http://rabbitvcs.org/">http://rabbitvcs.org/</a>

<sup>15 &</sup>lt;a href="http://pysvn.tigris.org/docs/WorkBench.html">http://pysvn.tigris.org/docs/WorkBench.html</a>




- 5. Download Eclipse IDE for GNU/Linux OS. As recommended in the Developer Handbook of universAAL, we will get the 3.5.2 version from Eclipse repositories<sup>16</sup>.
  - Extract the files on (for instance) ~/universAAL/eclipse/ (create the folder).
  - As a recommendation, create a new (clean) workspace folder. For example, ~/universAAL/uAAL\_workspace/. Each universAAL project will be a sub-folder of this one.
  - Start Eclipse IDE and select ~/universAAL/uAAL\_workspace/ as workspace for this session (this name is not important now, since we are only adding plugins to Eclipse but is a good practice).
    - 8. Install Eclipse required plugins:
  - Install Subversion plugin<sup>17</sup> for Eclipse IDE using the normal procedure.
    - Following the recommendation from the universAAL Developer Handbook, install a svn command line client. In terminal console with root privileges, enter the command sudo apt-get install subversion.
  - Install M2Eclipse plugin<sup>18</sup> for Eclipse IDE.

<sup>16</sup> http://archive.eclipse.org/eclipse/downloads/drops/R-3.5.2-201002111343/download.php?dropFile=eclipse-SDK-3.5.2-linux-gtk.tar.gz

<sup>17&</sup>lt;a href="http://subclipse.tigris.org/update\_1.6.x">http://subclipse.tigris.org/update\_1.6.x</a>

<sup>18</sup>http://m2eclipse.sonatype.org/sites/m2e

- Install Pax Runner plugin<sup>19</sup> for Eclipse IDE.
- 9. After installing all the additional software for Eclipse, the Installed Software window should look like this:




- 10. Configure Eclipse to look for the Java JDK.
  - Open the file eclipse.ini in ~/universAAL/eclipse/ folder.
  - Look for -wmargs line. Above it, write two lines (please keep in mind that the previous two lines have to be added as two lines, not one). Save and exit.

-vm /usr/lib/jvm/java-7-oracle/jre/bin/java


<sup>19</sup>http://www.ops4j.org/pax/eclipse/update/

- 11. Configure Maven options in the local system checking that the folder ~/.m2 is present.
  - Copy the settings.xml file available in **Resources** folder<sup>20</sup> (recently created at your computer) to your local .m2 folder.
- 12. Restart Eclipse. You may see a warning message saying that "Subversion 1.6 has a bug with GNOME keyring via the Subversion JavaHL API". Just press OK button.



13. Probably, you will have problems in order to have the JavaHL library working properly. This is used by Subversion as a high level API and serve as the JNI bridge between Java code and native libraries. You can check if the problem is present going to  $\mathbf{Window} \to \mathbf{Preferences} \to \mathbf{Team} \to \mathbf{SVN}$  (inside Eclipse IDE). If SVN Interface menu option shows "Not available", follow the next steps. If it shows JavaHL(JNI) 1.6.x, you should go without problems to the point number 14.

<sup>20~/</sup>universAAL/svn/support/trunk/resources/maven-repo-settings




 In GNU/Linux OS, the JavaHL library is in a separate package dependent on the main subversion package. To get it launch the command console and follow next steps:

```
sudo apt-get install libsvn-java
find / -name libsvnjavahl-1.so > ~/pathJavaHL
```


- Open the file eclipse.ini in ~/universAAL/eclipse/ folder.
- After -wmargs line, add the following line:
   -Djava.library.path=/usr/lib/jni/ (without quotes). This will create
   a file text with the path of the JavaHL library inside it. This information
   will be used in the next steps (typically, it will be /usr/lib/jni). Save
   and exit.

```
-startup
plugins/org.eclipse.equinox.launcher_1.0.201.R35x_v20090715.jar
--launcher.library
plugins/org.eclipse.equinox.launcher.gtk.linux.x86_1.0.200.v20090520
-showsplash
org.eclipse.platform
--launcher.XXMaxPermSize
256m
-vm
/usr/lib/jvm/java-7-oracle/jre/bin/java
-vmargs
-Djava.library.path=/usr/lib/jni/
-Xms40m
-Xmx256m
```

• Open again Eclipse IDE, and check if the SVN Interface option described before has a correct value.



14. At this point, you should be able to import any uAAL project as well as run it correctly.



#### **Troubleshooting**

1. The examples (and any universAAL bundle) will not run if the file sodapop.key is not in the proper folder. For solving this problem, look for any sodapop.key file in the ~/universAAL/svn/support/trunk/samples local folder, and copy it into [workspace\_folder]/rundir/mw.bus.model. If any of the folders in the later path does not exist, just create them.

#### ANDROID DEVELOPMENT SUPPORT

Setting up a real device for development/debbugging

http://developer.android.com/tools/device.html

```
angel@angel-laptop:~/Descargas/sw/Android/SDK_Google/android-sdk-linux/platform-tools$ sudo ./adb devices
[sudo] password for angel:
List of devices attached
emulator-5554 device
HT23XW136600 device
angel@angel-laptop:~/Descargas/sw/Android/SDK_Google/android-sdk-linux/platform-tools$
```

```
root@angel-laptop:/etc/udev/rules.d# ls -l
total 16
-rw-r--r-- 1 root root 71 sep 19 16:22 51-android.rules
-rw-r--r-- 1 root root 0 sep 19 16:21 51-android.rules~
-rw-r--r-- 1 root root 1055 sep 19 16:24 70-persistent-cd.rules
-rw-r--r-- 1 root root 629 abr 9 20:26 70-persistent-net.rules
-rw-r--r-- 1 root root 1157 oct 8 2011 README
root@angel-laptop:/etc/udev/rules.d#
```

